



SHOFU DIAMOND POINTS

Dental diamond burs

EN Instructions for use

CE 0123 04:2019-04
W000003

READ CAREFULLY BEFORE USE

INDICATIONS FOR USE

Finishing natural teeth and restorations

PRECAUTIONS

1. This product is intended for used by dental professionals only.
2. Do not use this product for any purpose other than specifically outlined in these instructions for use.

DIRECTIONS FOR USE

1. Sterilize diamond burs before each use.
 - Sterilization method: autoclave this product for 3 min at 134 °C (273 °F).

Cycle Type	Holding time	Temperature	Dry time
Pre-vacuum (3 times)	3'00"	134°C/ 273°F	10'00"

- Remove any contaminants, using washer disinfectant.
 - Do not immerse diamond burs into chlorine disinfectants.
 - Follow the instructions for use of washer disinfectant and sterilizer, respectively.
2. Mount a diamond bur on a dental handpiece (dental air turbine, micro motor etc.).
 - Insert the shank completely into the chuck of the handpiece, following the manufacturer's instructions.
 - Perform test run to check safety and smooth operation without run-out prior to use.
 - Make sure that head of a diamond bur is in good condition before use. Do not use damaged or defective instruments i.e. corrosion, surface flaw, deformation, dirt, exposure of shank due to repeated use.
 3. Use following the standard manner.
 - Never exceed the maximum permissible rotation speed.
 - Use local dust extractor, protective mask, etc while finishing to avoid adverse effect of dust to the human body.
 - Use protective equipment such as glasses while using.
 - Use intermittently with gentle touch and water spray.
 - Avoid any use at acute angles or with excessive pressure to thin and long abrasives. They are fragile and may bend or fracture under excessive loads.
 4. For reuse after using them intraorally, remove any contaminants immediately using washer disinfectant and sterilize.

PRODUCTS

- Robot Points
- Hybrid Points
- T&F Hybrid Points
- TPE Diamonds
- GingiCurettage



Manufacturer

SHOFU INC.

11 Kamitakamatsu-cho, Fukuine,
Higashiyama-ku, Kyoto 605-0983, Japan
www.shofu.com

EC REP Authorized Representative

SHOFU DENTAL GmbH

Am Brüll 17, 40878 Ratingen, Germany

SHOFU DENTAL CORPORATION

1225 Stone Drive, San Marcos, California
92078-4059, USA

**For the details about REPROCESSING
procedures, visit www.shofu.de**