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SHOFU Sep

Wax Separating Agent

SHOFU Sep is a wax separating agent intended for preventing wax from sticking to gypsum die, artificial teeth, metals, etc. Wax can be removed easily by applying this material before waxing up. It can be used with brush-on technique, dipping technique and drop-on technique.

DIRECTIONS FOR USE

1. Remove the polyethylene inside plug from the glass bottle. (The plug is for safe transportation, thus it is unnecessary for storage. The packing inside of the bottle cap should not be removed.)
2. Apply SHOFU Sep thinly to the separating surface with the brush provided.
Avoid applying excessively or insufficiently.
It is also possible to dip the model into SHOFU Sep, then excess material is removed.
3. Wax up onto the SHOFU Sep applied surface following the standard method.

NOTE ON USE

1. Thoroughly rinse waxed up reverse pin facing or pontic with water before cementation.
2. Apply fairly thinly before waxing up reverse pin facing, pontic, ceramic, metal or plastic for easy removal.
3. Exposure to direct sunlight might discolor the material. However, this does not affect the performance of this product.
4. Tightly close the cap after each dispensing.

INDICATIONS FOR USE

1. Separation of wax from dies of a variety of materials
2. Separation of wax from the cavity surface for direct wax pattern method
3. Separation of wax when taking wax patterns of pin holes or post crown holes of reverse pin facing or porcelain pontic, etc

PACKAGE

SHOFU Sep: 50 g (with a brush)

Manufacturer
SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan
SHOFU DENTAL TRADING (SHANGHAI) CO., LTD. No. 645 Jiye Road, Sheshan Industrial Park, Songjiang, Shanghai 201602, China
SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12, The Alpha, Science Park II, Singapore 117684