

ZIRCONOMER *P*

Zirconia Reinforced Restorative



ZIRCONOMER *P*

Zirconia Reinforced Restorative (ZRR)

Built on zirconia reinforced glass ionomer technology, **ZIRCONOMER *P*** is a self-adhesive, bulk-fill, nano-zirconia reinforced restorative material developed for strong and durable posterior restorations.

Functional Restorations for All Ages...

Class I & II
Restorations

Pediatric
& Geriatric
Restorations

Strong
Structural
Core
& Base

Silver Modified
Atraumatic
Restorative
Technique
(SMART
Technique)

Long-term
Transitional
Restorations

Atraumatic
Restorative
Treatment
(ART)

Biocompatible
Amalgam
Alternative





STRONG

ZIRCONOMER P exhibits a robust structure that strengthens and reinforces tooth structure while providing a secure marginal seal, and is considered an ideal biocompatible amalgam alternative for lasting posterior restorations. (compressive strength – more than 300 MPa)



SIMPLE

Reduce chair time – no priming, bonding or extensive tooth preparation!

Just mix, fill and finish at one go, to bulk fill cavities of all shapes and sizes while preserving maximum healthy tooth structure.



STABLE

No light curing required!

The self-adhesive restorative material eliminates the risk of polymerization shrinkage and displays a strong, stable bond to natural tooth structure with minimal interfacial stresses to assure greater longevity with minimal post-op sensitivity and optimal dimensional stability.



SAFE

Biocompatible with no risk of mercury!

A safe amalgam alternative, **ZIRCONOMER P** protects against secondary caries with its fluoride rich environment. Suitable for standard and high caries risk cases, the bulkfill, self-adhesive material restores a wide range of cavities in the posterior zone with ease and predictability.

SHOFU Smart Solutions



Prepare



Mix



Fill



Finish



Coat

1. Minimally invasive approach to cavity preparation

2. Mixing time
⌚ 30"
Mixing Ratio
Powder : Liquid
(2:1)

3. Working Time
(23 °C)
From start of mixing
⌚ 1'30"

Setting Time
(37 °C)
From end of mixing
⌚ 3'00"

4. Finish & Coat

Order Information



Zirconomer P 1-1 Set [REF C0040]
(Powder 12g & Liquid 5ml)

Related Product



Coat-It
[REF Y0055]



Composite Polishing Kit CA
[REF 0310]



SHOFU DENTAL ASIA-PACIFIC PTE. LTD.
1 Science Park Road, #01-01 / 02 The Capricorn, Singapore Science Park II,
Singapore 117528 Tel: (65) 6377 2722
eMail: mailbx@shofu.com.sg www.shofu.com.sg



ZIR20230930