





FOR MORE INFO

#### STREAMLINING ADHESIVE SOLUTIONS

**BeautiLink SA** is a two-component, dual-cure self-adhesive resin cement designed for swift, simple, and secure adhesion to a wide range of dental substrates.

Its novel multifunctional monomer formulation with unique silane coupling agent ensures strong, durable adhesion in both light and chemical cure modes, with virtually no post-op sensitivity. Radiopaque, easy to clean, and shade-stable, BeautiLink SA streamlines clinical workflows and enhances daily aesthetics.



# Available in:

#### Automix Syringe 5ML (Mixing tip x 10 pcs)

- Clear [REF C0203]
- Ivory [REF C0204]
- Opaque [REF C0205]

### **Handmix Syringe 9ML**

- Clear [REF C0200]
- Ivory [REF C0201]
- Opaque [REF C0202]

## **CLINICAL APPLICATIONS**

- Crown, Bridge, Inlay, Onlay, Veneer
- Adhesive Bridges (Maryland)

#### **FEATURES**

- Strong & lasting bond to multiple substrates
  - · Zirconia, Lithium Disilicate, Hybrid Ceramics
  - Metal, Indirect Composite
  - Enamel & Dentin
- **Low film thickness** (14  $\mu$ m) for invisible margins
- Superior masking capability for discolored abutments
- Easy, mess-free, clean-up with tack-cure removal
- ✓ Remarkable Shade stability (3 shades: Clear, Ivory, Opaque)
- ✓ Virtually no post-op sensitivity
- Synergy with **BeautiBond Xtreme** for nonretentive restorations
- Incorporates S-PRG filler technology for anticaries and anti-bacterial benefits

### **DIRECTION FOR USE**

- **1. Pre-treatment:** Clean abutment, sandblast if required, rinse & dry.
- 2. Dispensing: Dispense equally.
- **3. Application:** Apply to adhesive surface, seat restoration with proper pressure.
- **4. Tack Cure & Clean-up:** Tack cure for 1–2 sec, remove excess at margins
- 5. Final Cure: Light-cure each surface (LED 10 sec / Halogen 20 sec) or allow self/chemical cure.

\*\*Opaque Shade: Leave 2-3 min for self-cure & remove excess cement. Allow 5 min for chemical cure after light curing.

Working Time: 2 min I Setting Time: 3 min

# **RELATED PRODUCTS**

■ Beautibond Xtreme (REF Y2448)



Official Partner

