

ACRYLIC RESIN TEETH

Simply Natural and Economical
Cross-Linked / Abrasion Resistant



ACRYLIC RESIN TEETH – Mould Chart



UPPER

Mould Number

Length of Central Incisor : mm (minus collar)

Width of Central Incisor : mm

Width of 6 Anteriors (flat) : mm

Shades: A1, A2, A3, B2

LOWER

Mould Number

Length of Central Incisor : mm (minus collar)

Width of Central Incisor : mm

Width of 6 Anteriors (flat) : mm

Upper Anteriors • Bio Form

■ SQUARE

424

10.0
8.0
44.5



425

10.0
8.5
47.5



446

9.0
9.0
49.5



■ TAPER

513

9.0
7.5
41.5



514

10.0
8.0
43.5



515

10.5
8.5
46.5



516

11.0
9.0
49.5



517

11.5
9.5
51.5



534

10.5
8.0
44.5



535

11.0
8.5
47.5



543

8.5
7.5
42.0



544

9.5
8.0
44.5



546

10.5
9.0
49.0



Convenient wax-free bulk pack, economical and easy to use. Individual pontic types are stored in separate compartments for easy identification.



■ OVOID

615 10.5 8.5 48.0		624 10.5 8.0 44.5		625 11.0 8.5 48.0	
633 8.5 7.5 41.5		634 9.0 8.0 43.5		635 9.5 8.5 45.5	

Lower Anteriors • Bio Form

■ LOWER

23 8.5 31.5		25 9.5 34.5		26 10.0 36.0	
27 10.5 38.5		33 8.5 31.0		34 9.0 33.5	
35 9.5 34.5		36 10.0 36.0			



ACRYLIC RESIN TEETH – Mould Chart



UPPER

Mould Number
Length of first premolar : mm (minus collar)
Width of 4 Posteriors

Shades: A1, A2, A3, B2

LOWER

Mould Number
Length of first premolar : mm (minus collar)
Width of 4 Posteriors : mm

Upper / Lower Posteriors • Bio Form

M TYPE

M28

8.0
28.0



M30

9.0
30.0
9.0
33.0



M32

10.0
32.0
10.0
35.0



S TYPE

S28

6.0
28.0
6.0
31.0



S30

7.0
30.0
7.0
33.0



S32

8.0
32.0
8.0
35.0



Order Information



Acrylic Resin Teeth Anteriors
Set of 6's (24 sets in a box)



Acrylic Resin Teeth Posteriors
Set of 8's (16 sets in a box)



Acrylic Resin Teeth Anteriors
Bulk Pack x 120's



Acrylic Resin Teeth Posteriors
Bulk Pack x 160's



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605 - 0983, Japan (www.shofu.co.jp)

SHOFU DENTAL ASIA-PACIFIC PTE. LTD. 10 Science Park Road, #03-12 The Alpha, Science Park II, Singapore 117684 (www.shofu.com.sg)

SHOFU DENTAL TRADING (SHANGHAI) CO. LTD. No.645 Jiye Road, Sheshan Industry Park Songjiang Shanghai 201602, China (www.shofu.com.cn)

SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA (www.shofu.com)

SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany (www.shofu.de)

SHOFU UK. Riverside House, River Lawn Road, Tonbridge, Kent TN9 1EP, United Kingdom (www.shofu.co.uk)