

Glaslonomer FX ULTRA

CLICK / SCAN



FOR MORE INFO

UNIVERSAL ENHANCED GLASS IONOMER RESTORATIVE

When aesthetics and strength are combined for predictable performance

- **Ideal Translucency** for anterior restorations
- **Strength and Durability** for posterior restorations

Unique new formulation exhibits high initial reactivity and rapid maturation process provides remarkably stable, strong and fluoride-rich restorations for universal application



Available in 4 Shades – A2, A3, A3.5, B2

- **Glaslonomer FX ULTRA, 1-1 Set**
Powder 15g, Liquid 8ml (10g),
Measuring Scoop, Spatula, Cocoa Butter,
Mixing Pad, Universal Strips
- **Glaslonomer FX ULTRA, 1-1 Mini Set**
Powder 10g, Liquid 5.2ml

CLINICAL APPLICATIONS

- Paediatric & Geriatric restorations
- Class III, V and intermediate Class I & II restorations
- Root caries
- ART & minimally invasive treatment
- Tooth wear-erosion, abrasion & wedge shape defects
- Core build-up & in sandwich technique

FEATURES

- ✓ **Fast & convenient**
 - Self-adhesive & bulk fill
 - Easy to mix with good wettability
 - Ideal working time and a sharp-set reaction
- ✓ **Strong & durable**
- ✓ **Enhanced aesthetics with shade stability**
- ✓ **Radiopaque & moisture tolerant**
- ✓ **Low solubility and acid erosion**
- ✓ **Protect against secondary caries**

DIRECTIONS FOR USE

- Powder / Liquid ratio – 1 scoop : 1 drop
- Working time – 2 min 20 sec
(from start of mixing)
- Nett setting time – 2 min 30 sec
(from end of mixing)

RELATED PRODUCTS

- Coat It – Light cure high gloss resin coat
- Composite Finishing Kit



Official Partner

