

BEAUTIFIL II Enamel

CLICK / SCAN



FOR MORE INFO

UNIVERSAL BIOACTIVE NANO-HYBRID VALUE-BASED COMPOSITE RESIN

Designed for clinicians who focus on minimally invasive cosmetic restorations with composite artistry to mimic natural enamel ...

- **4 Age Defined Shades** with translucency and opalescent to mimic natural enamel
- **Easy to Adapt** to create fine surface details
- **Complimentary Shade Line-Up** for Beautifil II, Beautifil II LS & Injectable X

Formulated to provide 4 value based shades with varying levels of brightness to mimic natural enamel and create invisible restorations that offer predictable clinical performance and long-term durability.



Available in 4 Shades – 2.5gm syringe



CLINICAL APPLICATIONS

- Highly aesthetic direct veneers
- Post Ortho Incisal adjustments
- Cosmetic tooth modifications to enhance smiles
- Enamel layer for all types of restorations
- Root caries

FEATURES

- ✓ **Aesthetics that mimic natural Enamel**
 - Natural opalescence & translucency
 - Shade stability before and after cure
 - Easy to polish with sustained polish retention
- ✓ **Bioactive to protect against secondary caries**
 - Minimize plaque formation
 - Reduce cariogenic bacteria
- ✓ **Superior handling**
 - Non-sticky with smooth creamy consistency
 - One push syringe for controlled delivery
 - Long working time to easily create fine surface details
- ✓ **Lasting restorations**
 - Remarkable stain & wear resistance
 - Strong and durable

PROPERTIES

- Flexural Strength: 126 Mpa
- Filler Load: 82.6 wt%
- Radiopaque

RELATED PRODUCTS

- Beautifil II/ Beautifil II LS / Beautifil Injectable X
- BeautiBond / FL Bond II
- Super-Snap Xtreme Disks / OneGloss 2-in-1 Smart Polisher



Official Partner

