



FOR MORE INFO

# **UNIVERSAL BIOACTIVE** NANO-HYBRID COMPOSITE RESIN

Designed for clinicians who prefer a non-sticky, easy to adapt composite with a smooth consistency that offers ...

- Colour Stability & Easy Shade Match for life-like aesthetics
- Strong and Durable for predictable clinical performance
- Bioactive Properties for prevention of secondary caries

The patented bioactive S-PRG filler technology with a unique combination of six beneficial ions helps provide bioactive properties, while maintaining aesthetics, physical, mechanical and optical properties of direct resins.

Simplified shade line-up helps to minimize inventory while achieving the desired aesthetic outcomes.



#### Available in 14 Shades - 4.5gm syringe

- Universal Dentin A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3
- Opaque Dentin A2O, A3O
- Enamel INC, Bleach White BW

### **CLINICAL APPLICATIONS**

- Direct cosmetic tooth modifications, Diastema closure
- Class I VI cavities & root caries
- Tooth wear abrasion, erosion, attrition, abfraction
- Repair of fractured incisal edges
- Direct veneers and core build-up

### **FEATURES**

### True to nature aesthetics

- 3 different opacities with natural opalescence
- Blends well even with a single shade
- Shade stability before and after cure

# Protect against secondary caries

- Minimize plaque formation
- Neutralize the sugar acids
- Reduce cariogenic bacteria

## **Superior handling**

- Non-sticky, non-tacky and smooth consistency
- Easy to adapt and place

#### Long lasting restorations

- Remarkable wear resistance
- Clinically proven performance with 8 and 13 year clinical trials

## **PROPERTIES**

Flexural Strength: 130 Mpa

Wear Resistance: 0.52 wt%

Filler Load: 83 wt%

### **RELATED PRODUCTS**

- Beautifil II Enamel & Gingiva
- BeautiBond / FL Bond II
- Composite Polishing Kit
- Super-Snap Xtreme Disks / OneGloss 2-in-1 Smart Polisher



