

# CERAMAGE

CLICK / SCAN



FOR MORE INFO

## ZIRCONIUM SILICATE INDIRECT COMPOSITE SYSTEM

**Designed to meet patient-specific aesthetic and functional needs ...**

- Life-like aesthetics & colour reproduction
- Strong, Durable & Resilient with excellent shock-absorbing capability
- Long-term colour stability with minimal abrasion of the antagonist and easy intra-oral adjustment
- High wear and stain resistance

Ceramage light-cured indirect composite is formulated with a high filler load and a robust resin base for a range of removal and fixed prosthesis and routinely prescribed for implant-supported prosthesis.



**Available in kits and refills:**

- Ceramage Standard Set (REF Y1806U)
- Ceramage Gum Color Full Set (REF Y1804U)
- Refills – Paste in 4.6gm syringe / Flowable in 2 ml

## CLINICAL APPLICATIONS

### ■ Fixed:

- Metal-free / Reinforced Crowns
- Metal / Fiber Reinforced Bridges
- Implant-supported Crown & Bridges
- Inlays / Onlays / Veneers
- Long-term Temporary Restorations
- Repair of porcelain / Indirect Resin Restorations
- Modification of Artificial Resin Teeth

### ■ Removable:

- Overdentures
- Precision Attachment Work
- Removable Telescopic & Implant-supported Restorations
- Complete / Partial Dentures
- Individualized Gum in Removal Prosthesis

## FEATURES



### Natural Aesthetics

- Like-life restorations with natural opalescence & perfect shade match
- Individual characterization to replicate patient's unique dentition



### Remarkable Mechanical Properties

- Strong, flexible & durable with high colour stability



### Optimal Handling Characteristics

- Easy manipulation, non-sticky with ideal viscosity
- Paste & flowable variants



### Simple & Fast to Fabricate

- Effortless to build-up in precise layers
- Easy to polish to a natural gloss



### Shock-resistant & Repairable

- Resilient & Functional restorations

## PROPERTIES

- Compressive Strength: 354 Mpa
- Flexural Strength: 146 Mpa
- Filler Load: 73 wt%

## RELATED PRODUCTS

- Porcelain Primer / Universal Primer / HC Primer / CeraResin Bond Set
- Lite Art Resin Stains
- Ceramage Polishing Kit
- MiCD a-Cure LED

