

SHOFU BLOCK SHOFU DISK

CLICK / SCAN



FOR MORE INFO

DESIGNED TO PROVIDE A SEAMLESS RESTORATIVE SOLUTION FOR THE DIGITAL WORKFLOW WITH...

- Natural aesthetics
- Intra-oral modifications
- Flexibility & Shock resistance
- Easy processing and time-efficient
- Superior flexural strength

Developed with innovative filler technology and optimal filler-resin balance, SHOFU Block/ Disk HC hybrid ceramics offer life-like aesthetics, superior mechanical properties & faster milling for predictable restorations, with the option of same day delivery.



Available in:

- Sizes: Block – M & S
Disk – 98mm x 14mm
- 2 types of Pins – Universal & CEREC
- 11 Shades: HT (High Translucency) – 3 shades
LT (Low Translucency) – 6 shades
Enamel shade – 59, OC

CLINICAL APPLICATIONS

- Laminate Veneers
- Inlays / Onlays / Overlays
- Anterior & Posterior Full Crowns
- Implant Supported Restorations

FEATURES

- ✓ **Natural Aesthetics**
 - Natural opalescence & fluorescence
 - Remarkable polishability
 - Outstanding colour stability
 - Easy individual characterization
- ✓ **Enhanced Efficiency & Productivity**
 - Faster milling
 - Better accuracy & precision
 - No firing / sintering required
- ✓ **Predictable Clinical Performance**
 - Strong & durable
 - Greater marginal integrity
 - Resilience & flexibility
 - Greater wear resistance
 - Kinder to antagonist
 - Intra-oral adjustability & repairability

PROPERTIES

- Flexural Strength: 191 MPa
- Compressive Strength: 390 MPa
- Vickers Hardness: 66 Hv0.2

RELATED PRODUCTS

- Lite Art Resin Stains
- Resin Glaze Liquid & Resin Glaze Primer
- Ceramage Polishing Kit
- CompoMaster Coarse / CompoMaster
- Diamond Polishing Paste with Super-Snap Buff
- HC Primer & CeraResin Bond Set
- ResiCem
- Beauticem SA

