



## ALL PURPOSE TRUE-COLOUR STAIN SYSTEM

Designed for individual characterization of all types of ceramics to achieve life-like 2D or 3D stain effects...

- Low-fusing stains compatible with all conventional and CAD/CAM dental ceramics
- Ultra-fine glass technology assures easy, uniform application
- Intense colour pigment structure with natural fluorescence
- True-colour liquid to visualize results in real time

Formulated as an innovative, universal, vivid stain system that ensures predictable staining and characterization of all types of ceramics.



### Available in sets and refills:

- Vintage Art Universal Basic Set
- Stain Refills – 2gm
- Glazing Powder GP – 15gm, 50gm
- Glazing Powder GP-F – 15gm, 50gm
- Vintage Art Universal Liquid, 50ml
- Yamamoto Liquid, 50ml

## CLINICAL APPLICATIONS

- Staining of monolithic or layered ceramics
- Individualizing low & high fusing ceramics
- Ideal for both 2D & 3D staining
- Internal or external characterization

## FEATURES

- ✓ **Convenience of Use**
  - Microfine particle structure that assures ease of mixing and application in thin, even and bubble-free layers to create different contrasts
- ✓ **Universal application – Streamlined Inventory**
  - Compatible with monolithic and layered/ High or Low-fusing ceramics
- ✓ **Flexibility**
  - Option of 2 glaze powders – with or without fluorescence
- ✓ **Time-efficient, Predictable & Error-free**
  - True-colour (Yamamoto) liquid with glass-like refractive index to visualize shade effect and translucency of the stain mixture before firing

## FINAL FIRING TEMPERATURE GUIDELINE

- Veneered Lithium Disilicate  
730°C – 740°C
- Monolithic Lithium Disilicate  
740°C – 750°C
- Veneered Zirconia / PFM Restorations  
750°C – 800°C
- Monolithic Zirconia / Porcelain Teeth  
800°C – 850°C

## RELATED PRODUCTS

- Vintage Pro
- Vintage ZR
- Vintage LD
- SHOFU ZR Lucent
- SHOFU Porcelain Brush Set