



CLICK / SCAN



FOR MORE INFO

LITHIUM DISILICATE PORCELAIN

Designed to meet the highest level of aesthetic expectation with natural translucency and opalescence

- Predictable & flexible all-ceramic system with the option of press or layering ceramics
- Comprehensive shade line-up with exceptional colour match
- Remarkable strength & lasting durability
- Easy and quick divesting with minimal reaction layer on Vintage LD Press

Formulated with an innovative lithium disilicate based all-ceramic system that offers life-like aesthetics with consistent quality for long-term clinical performance.



Available in:

- Vintage LD Press Ingots – 3g each
4 Levels of Translucency – 30 shades
(T / MT / LO / MO)
- Vintage LD Porcelain sets and refills (15g, 50g)
Layering ceramics – 17 shades

CLINICAL APPLICATIONS

- **Restoration of Anterior / Posterior teeth...**
 - Veneers, Inlays, Onlays & Overlays
 - Anterior and Posterior Full Crown
 - 3 units Anterior Bridge
 - Implant Supported cases

FEATURES

- ✓ **Natural Aesthetics**
 - Efficiently simulates light transmission and diffusion of natural teeth for life-like results
 - Exceptional shade stability with impeccable colour match
- ✓ **Reliable & Flexible**
 - Simple fabrication – option of Press, Press & cut-back or Full build-up
 - Just 2 layers provides outstanding aesthetic outcome
 - Optimal precision for consistent results and fit
 - Quick, easy and clean divesting saves time (used with Ceravety Press & Cast)
- ✓ **Strong & Durable**
 - Strong, sturdy & integrated microstructure
 - Remarkable mechanical properties

PROPERTIES

- CTE: $10 \times 10^{-6} \text{ K}^{-1}$ (<25°C – 500°C)
- Glass Transition Temperature: 510°C
- Flexural Strength: 370 – 400 Mpa

RELATED PRODUCTS

- Ceravety Press & Cast
- Vintage Art Universal Stains
- SHOFU Pure Kolinsky Porcelain Brushes
- SHOFU Porcelain Brush Set
- ZII Master Finishing & Polishing Kit

