



CLICK / SCAN



FOR MORE INFO

THE AESTHETIC ELEMENT IN ALL-CERAMIC ZIRCONIA RESTORATIONS

Designed to create life-like all ceramic restorations with...

- Extensive range of Natural shades with opalescent & translucent enamel porcelain
- Feldspar-based ultrafine ceramics with dense, homogeneous structure for greater strength & durability
- Remarkable firing stability
- Suitable for all zirconia frameworks
- Excellent shade match

Micro-fine, high-fusing zirconia ceramics ideal for veneering all zirconia frameworks. Optimal light diffusion to achieve life-like shade harmonization to natural teeth even with thin layers.



Available in:

- Vintage ZR sets and refills (15gm, 50gm)
- 33 Shades: Standard Body – 16 shades
Incisal colours – 13 shades
Whitening colours – 4 shades
- Variety of Enamel Effect, Cervical, Cervical-Trans, Colour Effect and Gum Colours

CLINICAL APPLICATIONS

- Anterior & Posterior Full Crown
- Bridge and Implant supported restorations

FEATURES

- ✓ **Predictable Aesthetics**
 - Life-like aesthetics with a simple 2-layered build-up technique
 - Natural opalescence & fluorescence
 - Additional whitening colours
 - Corresponds to Vita shades
- ✓ **Enhanced Efficiency**
 - Simple & reliable layering of all standard zirconia frameworks
 - Easy to handle & build-up
 - Logical shade concept for easy application
- ✓ **Reliable Clinical Performance**
 - Homogeneous, dense surface with low firing shrinkage
 - Stable CTE, even with multiple firing
 - Greater wear resistance & less abrasion on opposing natural teeth
 - Excellent bonding characteristics with high fusing Opaque Liner

PROPERTIES

- CTE: $9.4 \times 10^{-6} \text{ K}^{-1}$ (25°C – 500°C)
- Solubility: 14.8 ($\mu\text{g} / \text{cm}^2$)

RELATED PRODUCTS

- Vintage Art Universal Stains
- SHOFU ZR Lucent Disk
- Porcelain Brush Set

