

Glaslonomer FX **ULTRA** CAPSULE

CLICK / SCAN



FOR MORE INFO

STRONG, PREDICTABLE & TRANSLUCENT

Highly translucent, self-adhesive and durable glass ionomer developed to meet the daily restorative needs of patients across the entire age spectrum, with the convenience of capsule delivery.

ADVANCED FX ULTRA GIC TECHNOLOGY

- Translucent & Aesthetic – Developed using highly translucent ‘HT’ Glass with dual surface treatment and anti-cluster molecules
- Strong & Durable - High-performance polycarboxylic acid co-polymer” fosters strong ionic bonds with multivalent ions to ensure long lasting restorations.
- One-step activation - deliver an easy, uniform & fast mix with optimal fluoride release & recharge



Available in:

Glaslonomer FX-Ultra Capsules, 50pcs
Shade : A2 (REF C0027) , A3 (REF C0028) ,
A3.5 (REF C0029)

CLINICAL APPLICATIONS

- Pediatric & Geriatric Restorations
Class III, V Restorations
- Erosions, Abrasions, Abfractions
- Core Build-up & Sandwich Technique
- Non-load bearing Class I, II Restorations
- Immediate Restorations in heavy stress bearing areas

FEATURES

- ✓ **Better Aesthetic Outcome**
 - Ultra-Translucent for easy shade match
- ✓ **Durable & Lasting**
 - Superior mechanical properties with greater wear resistance
- ✓ **Saves Chair Time**
 - Time-efficient with optimal working time & snap set reaction
- ✓ **Stable & Reliable**
 - Low solubility & acid erosion

SHOFU SMART SOLUTIONS

- 1. Prepare**
Minimally Invasive Cavity Preparation
- 2. Restore**
 - **Press to Active** (Do not click in capsule applicator)
 - **Mix** (10 sec mixing in capsule mixer)
 - **Fill** (Click 2 times and fill)
- 3. Finish & Protect**
Finish (Dura White stone) and apply COAT-IT

RELATED PRODUCTS

- Coat-it (REF Y0055)
- Composite Polishing Kit CA (REF 0310)

