

EyeSpecial C-V

CLICK / SCAN



FOR MORE INFO

SMART DENTAL CAMERA FOR DAILY PHOTOGRAPHY

The **EyeSpecial C-V** is a tailor-made dental camera designed to revolutionize clinical photography. Lightweight, compact, and intuitive, it captures consistent, true-colour images with minimal effort.

Its one-handed operation frees the other for dental tools, while the built-in autofocus and intelligent FlashMatic system deliver precise, glare-free photos in every shot. With hygienic, ethanol-resistant surfaces, it fully supports stringent cross-infection control protocols.



Available in:

E0028 EyeSpecial C-V (incl. accessories)

- E0022 Cross Polarizer
- E0002 Hand Strap
- E0004 Shofu Gray Card
- E0005 Lens Cap
- Case for Eye Special C-V

CLINICAL APPLICATIONS

- Daily dental photography in Clinical Practice for case documentation, patient/lab communication, medico-legal purposes
- Dental Lab Communication

FEATURES

- ✓ **9 Pre-set Shooting Modes** – Optimized for various dental applications (standard, low glare, whitening, tele-macro, isolate shade, movie, surgery, mirror, face).
- ✓ **FlashMatic System** – Patented auto-flash control for true-colour, glare-free photography
- ✓ **Continuous Autofocus** – Ensures sharp, consistent video recording.
- ✓ **One-Handed Compact Design** – Lightweight, ergonomic, and convenient for everyday use.
- ✓ **Intuitive Touchscreen** – Easy-to-navigate menus with icon-based guidance.
- ✓ **Enhanced Detail Capture** – Tele-macro mode with close-up lens support (up to 2:1 magnification).
- ✓ **Cross Polarizer (Optional)** – Eliminates enamel reflections to reveal internal tooth structures.
- ✓ **Durability** – Water- and chemical-resistant polymer body, easily disinfected with ethanol.

BRIEF SPECIFICATION

1. Focal Length – 50-217 mm (total 8 focal lengths)
2. Digital Zoom – 10 steps from x1.05 to x2
3. Maximum Aperture –
 - Wide: 4.6 – 5.3
 - Tele: 6.5 – 7.5
4. SD Card – SD/SDHC/SDXC; max.512GB

RELATED PRODUCTS

- Cross Polarizer (E0022)

