

# AI-INSTA

AI-Aided CAD/CAM

Same Day Digital Solution



Official Partner



# AI-Driven Prosthesis Within An Hour

Experience the Next-Generation of digital innovation with **AI-INSTA**. Powered by advanced AI (Artificial Intelligence), it features a fully automated, seamless CAD/CAM process that streamlines daily workflow from design to milling, delivering high quality prosthesis within an hour.

## No Complexity. No Overheads. Just Predictable Outcomes.

The innovative, AI-driven digital system ensures precision, reliability, efficiency and effortless adoption even for first time users, while saving time, labour and cost with pay per use subscription. **AI-INSTA** gives you total control with automated design and manufacturing of crowns, inlays, onlays, and veneers in the comfort of your clinic.



## Your All-in-One Digital Partner

**AI-INSTA** delivers a unified digital solution, designed to help clinics achieve consistent, high-quality results. From setup and training to service and support, **AI-INSTA** system performs at its best — every day.



**No software learning curve** — AI design and CAM processing are fully automated.



**No investment in high-end computer** — works without a dedicated laptop.



**No annual maintenance costs** — AI updates and new features applied automatically.



**On-demand cloud access** — through subscription for complete cost control.


## Smart Designing & Milling for an Efficient Chairside Digital Workflow

The AI technology provides an intuitive, step-by-step CAD/CAM function, making the AI-INSTA workflow simple and seamless.

All it takes is 3 mins for the AI CAD design, making it ideal for one visit treatment.

Just upload the patient's 3D scans and let the AI CAD design for you.

This simple to use, predictable cloud-based platform provides greater flexibility with hassle-free upgrading.



Select from an array of dental CAD/CAM materials to manufacture precision-engineered, accurate restorations that meet your patient's individual aesthetic and bio-functional needs.


### WET/DRY Milling of ...

- └ Hybrid Ceramics
- └ Glass Ceramics
- └ Zirconia
- └ PMMA



## The Ideal Material For Same Day Digital Restorations:

**SHOFU HC Hybrid Ceramic** is designed meticulously with a unique filler technology for exceptionally natural aesthetics, superior mechanical properties, faster milling and easy high-lustre polishing. Ideal for a spectrum of everyday aesthetic restorations, **SHOFU HC Hybrid Ceramic** can be easily milled with **AI-INSTA** to fabricate inlays, onlays, veneers, anterior/posterior crowns for same day delivery in your practice.




### Shade Range

	W2	A1	A2	A3	A3.5	B3	59	OC
<b>LT</b> Low Translucency								
<b>HT</b> High Translucency								
<b>Enamel</b>								


### Shape Dimensions (mm)

SHOFU BLOCK HC

Type	Dimension (mm)		
	a	b	c
<b>S</b>	10	12	16
<b>M</b>	12	14	18



SHOFU DISK HC



## Order Information



AI-INSTA [REF MST001]



### SHOFU BLOCK HC

5 blocks each per box; size S and M

### SHOFU BLOCK HC Intro Set

10 blocks, 1 each of sizes of S and M, shades A1-HT, A2-HT, A3-HT, A3-LT and OC

UNIVERSAL [REF 2150U]

CEREC [REF 2150C]



### SHOFU DISK HC

1 disk each per box; diameter 98mm, thickness 14mm


#### DISK HC 1-layer

SHADE	REF (1 pc.)
W2-LT	2180S
A1-HT	2186S
A1-LT	2181S
A2-HT	2187S
A2-LT	2182S
A3.5-LT	2184S
A3-HT	2188S
A3-LT	2183S
B3-LT	2185S
59 (Incisal)	2189S
OC (Occlusal)	2190S

#### BLOCK HC (UNIVERSAL)

1-layer

SIZE	SHADE	REF (5 pcs.)
M	W2-LT	2157U
	A1-HT	2158U
	A1-LT	2159U
	A2-HT	2160U
	A2-LT	2161U
	A3.5-LT	2162U
	A3-HT	2168U
	A3-LT	2169U
	B3-LT	2170U
	59 (Incisal)	2171U
OC (Occlusal)	2172U	
S	W2-LT	2151U
	A1-HT	2152U
	A1-LT	2153U
	A2-HT	2154U
	A2-LT	2155U
	A3.5-LT	2156U
	A3-HT	2163U
	A3-LT	2164U
	B3-LT	2165U
	59 (Incisal)	2166U
OC (Occlusal)	2167U	



## Related Products




**Lite Art Full Set [REF 1970]**

Contains: Colour Stains 1ml x 15, Clear Liquid 6ml x 1, Uni Brush (No.5) x 1, Colour Card x1, Nozzle Tip x 15, Nozzle Cap x 15



**Ceramage Polishing Kit HP**



**BeautiLink SA – Automix syringe 5ml**  
Mixing Tip X 10 Pcs

Clear [REF C0203]  
Ivory [REF C0204]  
Opaque [REF C0205]



**BeautiLink SA – Handmix syringe 9ml**

Clear [REF C0200]  
Ivory [REF C0201]  
Opaque [REF C0202]



**BeautiBond Xtreme Set [REF Y2448]**  
(5ml bottle, Micro brush 50's, V dish 25's)

Official Partner



**SHOFU DENTAL ASIA-PACIFIC PTE. LTD.**  
1 Science Park Road, #01-01 / 02 The Capricorn, Science Park II,  
Singapore 117528 Tel: (65) 6377 2722  
eMail: mailbx@shofu.com.sg www.shofu.com.sg



A120251130