

# ZirGloss

CLICK / SCAN



FOR MORE INFO

## ALL-CERAMIC SUPERPOLISHER FOR A GLAZED-LIKE SHINE

**ZirGloss** is a diamond-impregnated polishing paste designed for zirconia and lithium disilicate restorations.

It delivers a flawless, glazed-like surface finish, without the need for additional firing. Its dense, homogeneous blend of diamond and alumina particles in a beeswax base ensures excellent handling, controlled abrasive release, and minimal wastage.



Available in:

**ZirGloss [REF N0002]**  
20g polishing paste in a jar

### CLINICAL APPLICATIONS

- Super-polishing of all types of zirconia, lithium disilicate and porcelain restorations

### FEATURES

- ✓ **Easy** – Quick and effortless to super-polish to a high-gloss natural luster
- ✓ **Firm consistency** – splatter-free, easy to load onto Robinson brush/felt wheel
- ✓ **Economical** – minimal paste needed with low wastage
- ✓ **Residue-free** – easy wax removal from restorations
- ✓ **Low Antagonist wear** – Provides a satiny smooth, highly polished surface with very low roughness
- ✓ **Saves time** – Eliminates the need for re-firing after adjustments

### DIRECTION FOR USE

#### Predictable 3-Step Sequence:

1. **Contouring & Adjusting** – Dura-Green DIA
2. **Finishing & Pre-polishing** – ZilMaster Coarse/ Medium
3. **Final Polishing** – ZirGloss paste applied with Robinson's brush or felt wheel

### RELATED PRODUCTS

- **Dura Green DIA** (for contouring zirconia restorations)
- **ZilMaster Coarse & Medium** (for finishing and pre-polishing)

