

ZirGloss

Polishing compound for Zirconia,
Lithium Disilicate & Porcelain



ZirGloss

ZirGloss diamond impregnated polishing paste is designed to super polish **zirconia**, **lithium disilicate** and **porcelain** restorations. The smooth, densely packed, splatter-free, polishing paste is a perfect blend of diamond and aluminum oxide particles with ideal grit size, and a carrier paste of white beeswax.

The homogeneous and firm paste ensures controlled, pressure-responsive release of abrasives, to minimize wastage and ensure precise, high gloss polishing of all surfaces, greater restoration longevity, superior stain and abrasion resistance.



/ Time-Efficient

/ Economical, Easy & Efficient



/ Precise with Consistently Superior Polishing Results

Used in combination with Dura-Green DIA and ZilMaster, ZirGloss, guarantees a precise, lasting, high gloss surface in an easy, 3-step sequence.

01
Finish



Dura-Green DIA
Max. 30,000 rpm

02
Pre-Polish



For Zirconia-
ZilMaster Coarse
10,000-15,000 rpm
Max. 30,000 rpm



For Lithium Disilicate-
ZilMaster Medium
10,000-15,000 rpm
Max. 30,000 rpm

03
High
Gloss
Polish



ZirGloss



Order Information



ZirGloss [REF N0002]
(20g Polishing Paste in a jar)

Shape	CY5	FL3	IC5	KN7	RD3	TC2	WH6
ISO NO.	056	031	070	125	035	025	125
HP	0156	0158	0162	0154	0159	0161	0157

					
	Shape	Bullet	Minipoint	KN7	WH6
	ISO NO.	040	030	125	120
	ZilMaster Coarse HP 3 pcs	0640	0641	0642	0643
ZilMaster Medium HP 3 pcs	0646	0647	0648	0649	

Unmounted Wheels available



SHOFU DENTAL ASIA-PACIFIC PTE. LTD.
1 Science Park Road, #01-01 / 02 The Capricorn, Science Park II,
Singapore 117528 Tel: (65) 6377 2722
eMail: mailbx@shofu.com.sg www.shofu.com.sg



ZGL20251115